


Queen Mary, University of London
School of Physics and Astronomy

Useful Equations



- Equation-1: The conversion equation for energy (eV) and wavelength (nm):

$$E(\text{eV}) = \frac{1240}{\lambda(\text{nm})}$$

- Equation-2: The conversion equation for energy (eV) and temperature (K):

$$\frac{1\text{eV}}{k_B} = \frac{1.6022 \times 10^{-19} J}{1.3806 \times 10^{-23} J/K} = 11604.5052 K$$

- Equation-3: The X-ray absorption coefficient dependence on several parameters:

$$\mu \approx \frac{\rho Z^4}{AE^3}$$

where μ is the absorption coefficient, ρ is the sample density, Z is the atomic number, A is the atomic mass, and E is the X-ray energy.

- Equation-4: The pressure unit exchange relationship:

$$10\text{kbar} = 1\text{MPa}$$

$$1\text{MPa} = 10^6\text{Pa}$$

$$1\text{GPa} = 10^3\text{MPa} = 10^9\text{Pa}$$

- Equation-5: The commonly used formula to determine the lattice spacing in selected area electron diffraction (SAED):

$$R = \frac{L\lambda}{d}$$

The illustration for what each parameter means is given in the following figure: Here L and λ




Figure. 1 The illustration for SAED.

is constant, and therefore the measurement on the film plate R is directly linked to the lattice spacing d . The derivation is as following:

$$R = L \tan 2\theta \approx 2L \sin \theta$$



where θ is taken as quite small, thus we have $\tan 2\theta \approx 2\tan\theta \approx 2\sin\theta$. Then according to Bragg's law, we have:

$$2d\sin\theta = \lambda$$

By combining the above two equations we can obtain:

$$\frac{L}{d} = \frac{R}{\lambda} \Rightarrow R = \frac{L\lambda}{d}$$

This is just the commonly used formula for SAED pattern to extract lattice spacing information.

